

McConkey

EYE CARE

Scleral Lens Care

1

Soak the lenses in Boston Simplus over night
May also use Tangible Clean or ClearCare

2

When lenses are not "wetting" well, rub or store in Boston Simplus

3

Use saline to rinse lenses prior to inserting. Never use water

4

Use Scleralfil or other buffered saline to fill the lens before insertion.
May fill with Refresh Celluvisc to reduce midday fogging

5

Clean plungers with alcohol or soak in peroxide. Rinse before using

6

Use Refresh Relieva PF or any Optase drop for rewetting

Scleralfil is our preferred filling saline and vials should be replaced daily
Purliens bottles should be replaced every 10 days

Addipaks should be replaced daily, as they are preservative free.

Replace cases every 3 months (or more often)

Lenses without Hydrapeg coating may be stored dry, however, they must be sanitized each night

Please visit www.McConkeyEyeCare.com for additional resources, including videos on insertion and removal

Call or text: 443-453-5444
info@McConkeyEyeCare.com